


# HTE reaction planning

Before getting started:

- Select the file that is appropriate for your experimental design
  - Direct calibration curves and ee will require “ee\_noIS”
  - LCAP and >12 of something that isn’t preplated will require “lcap\_yield\_IS\_notpreplated”
  - Everything else (most things!) will require “lcap\_yield\_IS”
- Complete list of available ligands can be found by going to:  
<https://docs.google.com/spreadsheets/d/1AVvjK4NImj9yk8R7Be-ScvP-55ZgTwZerEuQbeGQft4/edit?usp=sharing>

# The underlying logic to the spreadsheet


- All experiments are described with sufficient detail to make troubleshooting and reproduction relatively easy
- All experiments are described such that translation of datasets into databases is relatively easy
- The basic version of this description is:



- Because inputting all of that information is challenging, the spreadsheet aims to streamline all of that


## 1. Adding Reaction Inputs: variable values

- All chemical and electrochemical inputs need to be added to the sheet, but that is split up between adding things to the plate on the “Plate Planning” sheet and on the “Reagents” tab
- Any input type where you are adding multiple variations should go on “Plate Planning”
- To enter the values, first select the role/role id for that column or row




# 1. Adding Reaction Inputs: variable values

- Once a role has been assigned, go ahead and add what the values are that you will be testing. This name can be as formal or informal as you'd like. If you are testing multiple concentrations or loadings and it's helpful to have that be explicitly in the label, then you should include it. These labels should be helpful to **you**



## 2. Adding Reaction Inputs: preplated ligands

- The ID numbers of any preplated materials can be found in the inventory linked on the first slide
- You can enter those ID numbers in the center of the plate as shown below



## 2. Adding Reaction Inputs: preplated ligands

- To confirm that the correct ligands have been added, you can check their names by going to the “Scan Ligands” tab
- It is recommended (but not required!) that you scan the barcodes on each tray as you are adding ligands, as confirmation that the correct one has been selected

The screenshot shows a spreadsheet application window titled "20240602\_Heck\_training". The spreadsheet has columns A through N. Columns A through G contain data for ligand trays. Column A is "well\_no", Column B is "ligand\_1\_no", Column C is "Barcode", Column D is "Number", and Column E is "Ligand Name". The data is organized into groups for trays A, B, C, D, and E. To the right of the main data table, there is a summary table titled "Ligand Tiers" with columns "Tier" and "Count".

well_no	ligand_1_no	Barcode	Number	Ligand Name
A1	461		461	DTBPF
A2	28		28	DavePhos
A3	33		33	DPPE
A4	34		34	DPPB
A5	3		3	P(tBu)3
A6	35		35	DPPF
B1	461		461	DTBPF
B2	28		28	DavePhos
B3	33		33	DPPE
B4	34		34	DPPB
B5	3		3	P(tBu)3
B6	35		35	DPPF
C1	461		461	DTBPF
C2	28		28	DavePhos
C3	33		33	DPPE
C4	34		34	DPPB
C5	3		3	P(tBu)3
C6	35		35	DPPF
D1	461		461	DTBPF
D2	28		28	DavePhos
D3	33		33	DPPE
D4	34		34	DPPB
D5	3		3	P(tBu)3
D6	35		35	DPPF

Ligand Tiers	Count
Tier 1	20
Tier 2	4
Tier 3	0
Tier 4	0
Tier 5	0
Tier 6	0
Tier 10	0
Tier 11	0

The bottom of the window shows a navigation bar with tabs: "Plate Planning", "Scan Ligands" (active), "Reagents", "Stock Solution Calculator", "Calibration Curves", "Uncorrected Area Counts", "Yields & LCAPs", and a "+" button. The status bar at the bottom indicates "Ready" and "Accessibility: Investigate".

### 3. Adding Reaction Inputs: electrochemical

- If you are testing multiple electrochemical conditions then can be entered in the same way
- If you are testing at one condition, then you should complete this section (24 well versions only)

product\_...

AutoSave ☐ 20240602\_Heck\_training

Home Insert Draw Page Layout Formulas Data Review View Automate

Comments Share

11 12 13 14 15 16

column variable base\_1

Preplated ligand ID numbers below.

row variable

row variable

47 NMeCy2

48 DIPEA

49 K3PO4

50 K2CO3

55

56

57

58

63

64

65

66

E-chem Conditions (leave blank if variable/not used)

Current (mA) 2

Voltage (V)

Charge (F/mol) 2

Anode Co

Cathode stainless

23 24 25 26 27 28

35 36 37 38 39 40

Dilution Factor

Reaction volume (μL):	100
Final reaction volume after dilution (μL):	300
First aliquot volume (μL):	25
First total analytical volume (μL):	500
Second aliquot volume (μL):	0
Second total analytical volume (μL):	0
Second aliquot volume (μL):	0
Second total analytical volume (μL):	0
Injection volume for calibration (μL):	2
Injection volume for data collection (μL):	2
Conversion factor (yield to final):	0.166667

General Notes

Research group:	
Experiment number:	
Analysis method:	
Stir rate:	
Time (h):	
Temperature (°C):	
Atmosphere:	Glovebox
Pressure (psi):	
Light source:	ambient
Light intensity:	
Solvent: (blank if variable)	
Internal standard:	TRUE
Reaction scale (μmol):	500
Preplated ligand loading:	0.2
Reaction vol (μL):	100
Ligand plating (μmol):	1

Tracked Compounds

Number	Name	SMILES	Yield (μmol)	Yield (mg)	[ ] (mg/mL)	Slope	R-squared
standard					0.000000	####	#DIV/0!
1					0.000000	####	#DIV/0!
2					0.000000	####	#DIV/0!
3					0.000000	####	#DIV/0!
4					0.000000	####	#DIV/0!
5					0.000000	####	#DIV/0!
6					0.000000	####	#DIV/0!
7					0.000000	####	#DIV/0!
8					0.000000	####	#DIV/0!
9					0.000000	####	#DIV/0!
10					0.000000	####	#DIV/0!

Correction Factors


Plate Planning Scan Ligands Reagents Stock Solution Calculator Calibration Curves Uncorrected Area Counts Yields & LCAPs

Ready Accessibility: Investigate 90%

- If you are not running an electrochemical reaction then leave these fields blank

## 4. Adding Reaction Inputs: non-chemical inputs

- There are many contributors to reaction to success beyond just chemical inputs and those are documented in the central column below the plate
- These inputs include stir rate, time, temperature, atmosphere, light source



- This is additionally where you can set the scale of the reaction with the variable that says “Reaction scale (μmol)”
- If you are using multiple main solvents then the Solvent field should be empty



## 5. Adding Reaction Inputs: everything else

- Once plate planning is complete, advance to the “Reagents” tab
- Any reagent that will be added to every well on a plate should be added to rows 1-10 as shown below
  - Concentration (M) should be used if the stock solution is not being prepared for the reaction in particular (i.e. a commercial solution, etc.)

	A	B	C	D	E	F	G	H	I	J	K	L	M
	non-variable	Code	Name	Source & Purity	Reactant Type	SMILES	FW (g/mol)	Loading (mol%)	Density (g/mL)	Concentration (M)	Loading per well (μmol)	Loading per well (mg)	Loading per well (p)
1	2	1	Bromobenzene	as found	substrate_1	BrC1=CC=CC=C1	157.01	100	1.5		20	3.1	2.0934
2	3	2	Methyl acrylate	as found	substrate_2	C=CC(=O)OC	86.09	200	0.95		40	3.4	3.624
3	4	3	Palladium acetate	Strem	catalyst_1	O=C(C)(C)C	224.51	10			2	0.4	
4	5	4									0		
5	6	5									0		
6	7	6									0		
7	8	7									0		
8	9	8									0		
9	10	9									0		
10	11	10									0		
11	column variable	11									0		
12	12										0		
13	13										0		
14	14										0		
15	15										0		
16	16										0		
17	column variable	23									0		
18	24										0		
19	25										0		
20	26										0		
21	27										0		
22	28										0		
23	column variable	35									0		
24	36										0		
25	37										0		
26	38										0		
27	39										0		
28	40										0		
29	row variable	47	NMeCy2		base_1						0		
30	48	DIPEA			base_1						0		

## 5. Adding Reaction Inputs: everything else

- Columns C-H are essential for every reagent including the variable ones entered on the previous sheet
- You will not be able to change any names of variable inputs here, but will need to switch back to “Plate Planning” to reassign those values
- Every reagent will need to have a unique label+id number. If they do not, then the cell will show an error

Code	Name	Source & Purity	Reactant Type	SMILES
1	Bromobenzene	as found	substrate_1	<chem>BrC1=CC=CC=C1</chem>
2	Methyl acrylate	as found	substrate_1	<chem>C=CC(=O)OC</chem>
3	Palladium acetate	Strem	catalyst_1	<chem>CC(=O)OC([Pd])(OC(=O)C)OC(=O)C</chem>

## 6. Figuring out stock solutions

- Most reagents will be added to the plate using stock solutions in the reaction solvent – this is the most straightforward and least error prone way of doing plate setup
- This means that for each well, the sum of the stock solution volumes should be equal to the desired reaction volume
- In general, you should aim to minimize the concentrations of stock solutions as much as possible as that will increase the accuracy
- In other words:
  - for a reaction run at 20  $\mu\text{mol}$  at 0.1 M, the total volume will be 200  $\mu\text{L}$
  - In this particular case, the stock solutions that need adding are:
    - $\text{Pd}(\text{OAc})_2$  – 2  $\mu\text{mol}$
    - Bromobenzene – 20  $\mu\text{mol}$
    - Methyl acrylate – 40  $\mu\text{mol}$
    - Base – 40  $\mu\text{mol}$
  - The 200  $\mu\text{L}$  can be divided evenly but if you have more than 4 stock solutions, then it is preferable to make the volumes of the higher loading components larger (i.e. more dilute)

## 6. Figuring out stock solutions

- To actually plan this out, go to the “Stock Solution Calculator” and enter the solvent, the number of wells the stock solution needs to go for, the desired stock loading, and the code for each reagent
  - The code is listed next to each name on “Reagents”
- This will automatically calculate how much you need to weigh out, the resulting concentration, and the amount of solvent you need to add
- If you need to do a slurry loading of a particular component then it is recommended to increase the number of wells to around 1.25x as many as are actually needed
- You do not need to be precise when weighing things out – if you list the actual mass weighed out then the sheet will account for that in the volume of solvent you need to add

Solution Number	Solvent	# of wells	Desired stock loading (µL)	First component Code	Name	Actual Mass (mg)	Amount for Stock (mg)	Amount for Stock (µL)	Stock Concentration (M)	Volume for Stock (mL)
1	DMA	24	50	1	Bromobenzene		94.2	62.8	0.4	1.437
2	DMA	24	50	2	Methyl acrylate		103.3	108.7	0.8	1.391
3	DMA	24	50	3	Palladium acetate		13.5	-	0.04	1.487
4					-		-	-		
5	DMA	24	50	47	NMeCy2		234.4	257.0	0.8	1.243
6	DMA	24	50	48	DIPEA		155.1	209.0	0.8	1.291
7	DMA	30	50	49	K3PO4		318.4	-	0.8	1.557
8	DMA	30	50	50	K2CO3		207.3	-	0.8	1.668
9					-		-	-		
10					-		-	-		
11					-		-	-		
12					-		-	-		
13					-		-	-		
14					-		-	-		

## 7. Running the experiment

- At this point you have entered all the plate inputs!
  - Hooray!!
- If you are on the 96-well version of the spreadsheet then there is a sheet that will color code wells according to how much of each reagent needs to be added

## 8. Analyzing the data

- Once you have finished running the experiment then the workbook will help you additionally figure out the dilutions required for your analytics and calculate yields
- To get started, go back to “Plate Planning” and complete the “Tracked Compounds” table with the names, SMILES, and yield in  $\mu\text{mol}$  of any compound present in your reaction mixture that you want to quantify
- If you have a correction factor for that reaction component then you can add it here as well

The screenshot shows the 'Plate Planning' tab in a spreadsheet application. The main area contains a grid for 'Reaction Conditions' with columns for 'Prepared Serial C numbers' and 'Yield (μmol)'. Below this, there are sections for 'Reaction Conditions (Serial Blank if variable not used)', 'Reaction Factor', 'General Notes', and 'Tracked Compounds'.

**Tracked Compounds Table:**

Number	Name	SMILES	Yield (μmol)	Yield (mg)	[ ] (mg/mL)	Slope	R-squared
standard	dodecane	CCCCCCC	3		0.000000		
1	Methyl cinnamate	CC(=C)C(=O)OC	10		0.000000	####	#DIV/0!
2	Bromobenzene	Brc1ccccc1	10		0.000000	####	#DIV/0!
3	Methyl acrylate	CC(=C)OC	20		0.000000	####	#DIV/0!
4					0.000000	####	#DIV/0!
5					0.000000	####	#DIV/0!
6					0.000000	####	#DIV/0!
7					0.000000	####	#DIV/0!
8					0.000000	####	#DIV/0!
9					0.000000	####	#DIV/0!
10					0.000000	####	#DIV/0!

**Correction Factors Table:**

Number	Name	SMILES	Yield (μmol)	Yield (mg)	[ ] (mg/mL)	Slope	R-squared
standard	dodecane	CCCCCCC	3		0.000000		
1	Methyl cinnamate	CC(=C)C(=O)OC	10		0.000000	####	#DIV/0!
2	Bromobenzene	Brc1ccccc1	10		0.000000	####	#DIV/0!
3	Methyl acrylate	CC(=C)OC	20		0.000000	####	#DIV/0!
4					0.000000	####	#DIV/0!
5					0.000000	####	#DIV/0!
6					0.000000	####	#DIV/0!
7					0.000000	####	#DIV/0!
8					0.000000	####	#DIV/0!
9					0.000000	####	#DIV/0!
10					0.000000	####	#DIV/0!

## 8. Analyzing the data

- Area counts can then be pasted into the “Uncorrected Area Counts” sheet but be sure to make sure that they are in the column with the heating correspond to the ID number listed on the “Plate Planning Sheet”
- In the 96-well version of the spreadsheet you will also see a notes column – if any well was not run as intended or not run, then this should be indicated here

AutoSave 20240221\_initialHeck Search (Cmd + Ctrl + U)

Home Insert Draw Page Layout Formulas Data Review View Automate

H5 543.9713

	A	B	C	D	H	L	P	T	X	AB
1	notes	well_no	total_areacount	standard	product_1_areacount	product_2_areacount	product_3_areacount	product_4_areacount	product_5_areacount	product_6_areacount
4	Run as intended	A1	247.026	4.3872						
5	Run as intended	A2	230.335	543.9713						
6	Run as intended	A3	231.952	337.1337						
7	Run as intended	A4	231.373	23.1081						
8	Run as intended	A5	225.627	7.8155						
9	Run as intended	A6	230.365	7.0371						
10	Run as intended	A7	229.908	200.4033						
11	Run as intended	A8	226.24	4.5189						
12	Run as intended	A9	223.618	116.2965						
13	Run as intended	A10	217.575	22.3279						
14	Run as intended	A11	225.085	46.9062						
15	Run as intended	A12	229.862	100.6428						
16	Run as intended	B1	223.151	13.4447						
17	Run as intended	B2	231.912	545.1774						
18	Run as intended	B3	226.073	153.7798						
19	Run as intended	B4	236.348	29.0075						
20	Run as intended	B5	228.645	11.827						
21	Run as intended	B6	214.866	1.4963						
22	Run as intended	B7	213.399	146.1039						
23	Run as intended	B8	215.558							
24	Run as intended	B9	218.062	88.698						
25	Run as intended	B10	220.32	24.4979						
26	Run as intended	B11	226.403	44.1985						
27	Run as intended	B12	206.737	82.4843						
28	Run as intended	C1	232.787	9.8682						
29	Run as intended	C2	212.214	481.5898						
30	Run as intended	C3	201.463	293.275						
31	Run as intended	C4	211.349	16.1824						
32	Run as intended	C5	221.03	5.593						
33	Run as intended	C6	207.674	29.1613						
34	Run as intended	C7	209.308	272.3993						
35	Run as intended	C8	217.622	298.5128						

Ready Accessibility: Investigate

Scan Ligands Reagents Stock Solution Calculator Calibration Curves Dosing Diagram **Uncorrected Area Counts** Yield +

100%



## 8. Analyzing the data

- If there is a correction factor and all yields are completed in “Plate Planning” then you should now see your yields in the “Yields and LCAP” tab

AutoSave

20240221\_initialHeck

Search (Cmd + Ctrl + U)

Home

Insert

Draw

Page Layout

Formulas

Data

Review

View

Automate

Comments

Share

D23

=VLOOKUP(A23,'Uncorrected Area Counts'!\$B\$4:\$AV\$99,7,FALSE)/VLOOKUP(A23,'Uncorrected Area Counts'!\$B\$4:\$AV\$99,3,FALSE)/VLOOKUP(1,'Plate'

	A	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q
1	well_no	one_lcap	one_yield	two_lcap	two_yield	three_lcap	three_yield	four_lcap	four_yield	five_lcap	five_yield	six_lcap	six_yield	seven_lcap	seven_yield	eig
4	A1	-1.77600738	0.7	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
5	A2	-236.165384	95.0	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
6	A3	-145.346134	58.5	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
7	A4	-9.98736242	4.0	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
8	A5	-3.46389974	1.4	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
9	A6	-3.0547596	1.2	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
10	A7	-87.1668121	35.1	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
11	A8	-1.99738862	0.8	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
12	A9	-52.0067705	20.9	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
13	A10	-10.2621672	4.1	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
14	A11	-20.8393274	8.4	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
15	A12	-43.7840857	17.6	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
16	B1	-6.02492839	2.4	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
17	B2	-235.079832	94.5	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
18	B3	-68.0221875	27.4	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
19	B4	-12.273221	4.9	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
20	B5	-5.17265887	2.1	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
21	B6	-0.69638752	0.3	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
22	B7	-68.4650645	27.5	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
23	B8	0	0.0	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
24	B9	-40.6756063	16.4	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
25	B10	-11.1192256	4.5	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
26	B11	-19.5220901	7.9	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
27	B12	-39.898257	16.0	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
28	C1	-4.23915607	1.7	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
29	C2	-226.936137	91.3	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
30	C3	-145.572853	58.5	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
31	C4	-7.65672308	3.1	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
32	C5	-2.53042573	1.0	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
33	C6	-14.0418366	5.6	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
34	C7	-130.14268	52.3	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	
35	C8	-137.170323	55.2	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	0.0	#DIV/0!	

Reagents

Stock Solution Calculator

Calibration Curves

Dosing Diagram

Uncorrected Area Counts

Yields & LCAPs

Ready

Accessibility: Investigate

</



## 9. Finishing up

- At this point your data can be visualized in notebook using the “Heatmap” tab, but you should also submit it to our visualization program to generate more complicated visualizations

